

ISO 489 Refractive Index Measurement System for Plastics â€™ Testing Equipment

Description

Designed for Refractive Index Determination According to ISO 489:2022

The **Refractive Index Measurement System** is designed for the characterization of plastic and polymeric materials based on the refractometric principles described in **ISO 489:2022 â€™ Plastics â€™ Determination of refractive index**.

The system is intended for applications in polymer laboratories, quality control, material identification, product development, and optical characterization of plastics.

ISO 489:2022 describes two measurement methods:

Method A â€™ Refractometric Method

Measurement using an Abbe refractometer or another refractometer capable of providing equivalent results.

Method B â€™ Immersion Method

Measurement using a microscope and immersion liquids of known refractive index.

For direct measurement of solid plastic specimens, sheets and films, the instrument is primarily based on the principles of **Method A**.

Table 2 — Immersion liquids

Immersion liquid	Refractive index at 23 °C n_D^{23}
n-Butyl carbonate	1,410
Tri-n-butyl citrate	1,444
n-Butyl phthalate	1,491
1-Bromonaphthalene	1,657
Diiodomethane (methylene iodide)	1,747
Aqueous solution of potassium mercury(II) iodide	1,419 to 1,733 ^a
Silicone oils	1,37 to 1,56 ^a
^a Useful range for the purpose of the test.	

How Does the Measurement Work?

The measurement is based on the optical interaction between the **plastic specimen** and a **high-refractive-index measuring prism**.

A small amount of contacting liquid is placed between the specimen and the main prism to establish proper optical contact.

The optical system illuminates the prism/specimen interface and generates a transition between a **bright region and a dark region**.

The angular position of this boundary depends on the critical optical condition at the interface and therefore on the **refractive index of the specimen**.




Figure 1 — Method for measuring refractive index of transparent sheet

The illustration should show:

- Light source
- Main prism
- Sub-prism
- Plastic specimen
- Contacting liquid
- Optical path
- Reflected/refracted light
- Detection system

Automatic Optical Boundary Detection

In a conventional Abbe refractometer, the operator observes the optical field through an eyepiece and manually aligns the boundary between the bright and dark regions with a reference crosshair.

A digital measurement system can replace this visual operation with an **electronic optical detector or imaging sensor**.




Figure 2 — Refractometer field of vision

The detector captures the intensity distribution across the optical field. The measurement algorithm determines the precise location of the light/dark transition and converts its position into a refractive-index value using the instrument calibration function.

This architecture can reduce operator dependency and enables automatic digital measurement.

The principle corresponds to the refractometer field-of-view arrangement illustrated in Figure 2 of ISO 489:2022.

Measurement Range

The Abbe refractometer described for Method A covers approximately:

Refractive Index

$n = 1.300$ to 1.700

Accuracy

$\hat{A} \pm 0.001$

This range covers many common polymeric materials and makes the system suitable for laboratory and industrial plastics characterization.

Optical Wavelength

Refractive index depends on the wavelength at which the measurement is performed.

ISO 489:2022 refers to the sodium D-line at approximately:

$$\lambda = 589 \text{ nm}$$

as the reference wavelength.

A modern digital instrument can therefore employ a controlled narrow-band optical source operating near **589 nm** to establish a defined measurement wavelength.

This is particularly important when comparing results between different instruments or laboratories.

Precision Temperature Control

Temperature has a direct influence on refractive-index measurement.

For Method A, determination is performed at:

$$23.0 \pm 0.5 \text{ }^{\circ}\text{C}$$

Therefore, temperature stability of the specimen and measuring prism is a critical part of the instrument design.

Recommended Temperature-Control Architecture

Temperature Sensor -> Digital Controller -> Heating/Cooling Element -> Measuring Prism

A digital instrument may use:

- Precision temperature sensor
- Closed-loop temperature controller
- Thermoelectric heating/cooling element
- Thermal coupling to the main prism
- Real-time temperature monitoring

The measurement should only be accepted when the prism/specimen assembly has reached the required temperature stability.

Measurement of Transparent Plastic Sheets

For transparent sheet specimens, the sample is placed directly against the measuring prism using a suitable contacting liquid.

The contacting liquid establishes an optically continuous interface between the plastic specimen and the prism.




Figure 3 — Method for measuring refractive index of film

Recommended cross-sectional arrangement:

Plastic Specimen

at

Contacting Liquid

at

Main Measuring Prism

The complete optical diagram should additionally show the illumination source, sub-prism and optical detection path.

The specimen surface contacting the prism should be optically flat and sufficiently well polished to provide reliable optical coupling.

For sheet specimens, ISO 489:2022 recommends a specimen thickness of approximately:

3 to 5 mm

Surface defects, cutting burrs, contamination and poor optical contact can adversely affect the measurement.

Measurement of Thin Plastic Films

The system can also be configured for measurement of very thin polymer films.

For film measurement, a small amount of contacting liquid is first placed on the main prism.

The film specimen is then positioned over the prism.

A glass plate is placed above the film using an additional layer of contacting liquid.




Figure 4 — Film with multiple refractive indices

Recommended cross-section:

Glass Plate

“

Contacting Liquid

“

Polymer Film

“

Contacting Liquid

at“

Main Prism

The arrangement provides optical coupling on both surfaces of the film and facilitates measurement of thin specimens.

ISO 489:2022 indicates that film thickness can be as low as approximately:

2 ¼m

Air bubbles between the film, liquid, glass and prism should be avoided because they can disturb the optical boundary and reduce measurement reliability.

Contacting Liquid

Selection of the contacting liquid is an important part of the measurement.

The liquid should have a **higher refractive index than the plastic material being measured.**

At the same time, it must not:

- Soften the polymer
- Dissolve the specimen
- Chemically attack the specimen
- Produce an undesirable interaction with the plastic

The liquid provides the optical connection between the specimen and measuring prism.

The appropriate liquid therefore depends on both the expected refractive index and the chemical compatibility of the polymer.

Table 1 — Contacting liquids

Plastic material	Contacting liquid
Cellulose derivatives	Aniseed oil or 1-bromonaphthalene
Fluorine-containing polymers	1-Bromonaphthalene
Urea-formaldehyde	Aniseed oil or 1-bromonaphthalene
Phenol-formaldehyde	1-Bromonaphthalene
Polyethylenes	1-Bromonaphthalene
Polyamides	1-Bromonaphthalene
Unsaturated polyester	1-Bromonaphthalene
Polyisobutylene	Saturated aqueous solution of zinc chloride made slightly acid
Poly(methyl methacrylate)	Saturated aqueous solution of zinc chloride made slightly acid or 1-bromonaphthalene
Polystyrene	Saturated potassium mercury(II) iodide solution
Styrene-acrylonitrile copolymers	1-Bromonaphthalene
Vinyl resins (vinyl chloride copolymer or plasticized PVC)	1-Bromonaphthalene
Poly(vinyl chloride)	1-Bromonaphthalene
Poly(ethylene terephthalate)	Methylene iodide
Polycarbonate	Methylene iodide
Diethylene glycol bis(allyl carbonate) (CR 39)	Methyl salicylate, aniseed oil or 1-bromonaphthalene
Polyarylate	Saturated aqueous solution of zinc chloride made slightly acid, methylene iodide or 1-bromonaphthalene
Polyetheretherketone	Methylene iodide
Polypropylene	1-Bromonaphthalene

Measurement of Anisotropic Plastics

Many polymer manufacturing processes create molecular orientation within the material.

Examples include:

- Extrusion
- Injection moulding
- Film stretching
- Oriented sheet production

As a result, the refractive index measured in one direction can differ from the refractive index measured in another direction.

ISO 489:2022 therefore includes procedures for measuring anisotropic materials in different specimen orientations.




Measurement direction	Thickness (Z axis)	At right angle to machine direction (X axis)	Thickness (Z axis)	Machine direction (Y axis)
Positioning of specimen with multiple refractive indices (machine direction indicated by arrows)	Light 	Light 	Light 	Light 
Polarizing direction (indicated by arrows)				

Figure 5 — Combination of specimen-positioning direction and polarizing-filter direction

The illustration should show a rectangular polymer specimen with:

X — Direction perpendicular to machine direction

Y — Machine direction

Z — Thickness direction

This diagram should be based conceptually on the orientation system shown in Figure 4 of ISO 489:2022.

Polarized Measurement

For anisotropic materials, the optical system can additionally incorporate a **polarizing filter**.

Measurements can then be performed using different combinations of:

- Specimen orientation
- Machine direction
- Transverse direction
- Thickness direction
- Polarization direction

Recommended illustration:

589 nm Light -> Polarizer -> Prism -> Oriented Polymer Specimen -> Detector

The polarization direction should be clearly indicated with arrows.

Such measurements are useful for investigating direction-dependent optical properties produced by polymer processing and molecular orientation.

Measurement of Translucent, Coloured and Opaque Plastics

The ISO 489 procedure is not limited to completely transparent specimens.

Certain:

Translucent “ Coloured “ Opaque

plastic materials can also be evaluated.

For these specimens, a **reflection measurement configuration** is used.

Instead of relying on transmitted light through the complete specimen, the optical system evaluates reflected light from the specimen/prism interface.




Figure 6 — Method for measuring refractive index of translucent, coloured and opaque materials

Recommended optical arrangement:

Light Source -> Illumination Window -> Main Prism -> Contacting Liquid -> Opaque Specimen

followed by:

Reflected Light -> Main Prism -> Optical Detection System

The diagram should clearly distinguish the **incident light** from the **reflected measurement light**.

This arrangement is conceptually based on the reflection configuration illustrated in Figure 6 of ISO 489:2022.

Sample Preparation

Correct specimen preparation is essential for accurate refractive-index measurement.

For sheet specimens, the surface in contact with the prism should be:

- Optically flat
- Sufficiently polished
- Clean
- Free from cutting burrs
- Free from contamination

The specimen should provide uniform contact with the measuring prism.

Poor surface preparation can produce an unstable or poorly defined optical boundary.

Environmental Conditions

Specimen conditioning is performed according to the conditions referenced in ISO 489:2022 and ISO 291.

The reference atmosphere specified is:

Temperature

$23 \pm 2 \text{ }^{\circ}\text{C}$

Relative Humidity

$50 \pm 5 \text{ \% RH}$

The measurement equipment should also be operated in a controlled environment to minimize variations caused by external temperature and humidity changes.

Digital Instrument Architecture

A modern refractive-index measurement system based on the ISO 489 Method A principle can integrate the complete measurement process into a digitally controlled platform.




Figure 7 — Becke line and its movement

Recommended system diagram:

589 nm Light Source

â†“

Illumination Optics

â†“

Sub-Prism / Sample Interface

â†“

Plastic Specimen

â†“

Contacting Liquid

â†“

Main Measuring Prism

±

Optical Detection System

±

Image / Signal Processing

±

Boundary Detection Algorithm

±

Calibration Function

±

REFRACTIVE INDEX (n)

A separate closed-loop temperature-control system should simultaneously maintain the measuring assembly at the required temperature.

Recommended System Features

A digital refractive-index measurement instrument can include:

- Measurement range approximately **1.300 to 1.700**
- Target measurement accuracy of **±0.001**
- Optical source near **589 nm**
- Precision optical measuring prism
- Digital optical detection
- Automatic light/dark boundary detection
- Closed-loop prism temperature control
- Measurement at **23.0 ± 0.5 °C**
- Digital calibration
- Sample identification
- Measurement history
- Test-data storage
- Automatic test-report generation
- USB or network data transfer
- Optional polarization module
- Optional reflection measurement mode

Typical Applications

The system can be used for optical characterization and quality control of:

- Transparent plastic sheets
- Polymer films
- Injection-moulded components
- Extruded plastics
- Oriented polymer films
- Translucent plastics
- Coloured polymers
- Selected opaque plastics

Typical application areas include:

Polymer Laboratories | Plastics Manufacturing | Quality Control | R&D | Material Identification | Incoming Material Inspection

Test Report

For traceable testing according to ISO 489:2022, the measurement report should include relevant information such as:

- Reference to ISO 489:2022
- Complete identification of the tested material
- Measurement method used
- Type of light source
- Measurement wavelength
- Specimen orientation where applicable
- Dispersion information where applicable to Method A
- Immersion liquid where Method B is used
- Measured refractive index
- Any deviation from the specified procedure
- Date of measurement

This provides a structured basis for laboratory documentation, material comparison and quality-control records.

Category

1. Equipment for Standards
2. Standards