



## Laboratory Injection Molding Machine-155mm

### Description

- Proper for different plastic materials //
- HMI 3â€³ touchscreen //
- PLC controlled //
- Cavities and sample shapes as per customer request will be quoted //
- The process is fully automatically controlled //
- Easy assembly of cavities //
- Suitable size of distance between tie bars for placement of molds for samples of Tensile, HDT, Vicat, Charpy, Izod, â€¦. //
- Molds for samples will be quoted separately based on customer request //
- Distance between tie bars 155mm //
- Cooling type is by air fan